



## Bleach Response Red Wine Stain

Studies were undertaken at 30°C, 40°C and 60°C to determine to bleach response for fresh stains supplied by Warwick Equest.


Results shown below are for 5cm stain circles, prepared on 100% knitted cotton, white fabric.

Wash Conditions:

- European Machine – Miele W487
- Normal wash, cotton cycle
- 3.5Kg clean ballast load
- 10g/l detergent dose
- Formulation contained 85% IECA\*, and where required 10% PCS and 5% Activator

Warwick Equest Red Wine stains show bleach response at 30°C, 40°C and 60°C compared to base detergent alone formulations. Red Wine stains are affected by pH in the absence of bleach systems and this is evident by a decrease in stain removal when using a base detergent (IECA\*) compared to water alone wash.

### Response for Red Wine



**Red Wine stains available include:** Red Wine, Bordeaux, Cabernet Sauvignon

For more details relating to stain selection and options, please contact us:

[enquiries@warwickequest.com](mailto:enquiries@warwickequest.com)

Ref: EQUEST BLEACHABLE #01.1 – RED WINE